

## GREENS DHS BARREL PLATE WITH DCP HOLE 135°- SHORT DIA BARREL [PEDIATRIC PLATES]



CAT. NO.SS	GREENS DHS Barrel Plate with DCP Hole 135°- Short Dia Barrel [Pediatric Plates]
288.102	Greens- DHS Barrel Plate with DCP Hole 135° Short Dia Barrel - 02 Holes, L- 46mm
288.103	Greens- DHS Barrel Plate with DCP Hole 135° Short Dia Barrel - 03 Holes, L- 62mm
288.104	Greens- DHS Barrel Plate with DCP Hole 135° Short Dia Barrel - 04 Holes, L- 78mm
288.105	Greens- DHS Barrel Plate with DCP Hole 135° Short Dia Barrel - 05 Holes, L- 94mm
288.106	Greens- DHS Barrel Plate with DCP Hole 135° Short Dia Barrel - 06 Holes, L- 110mm
288.107	GREENS DHS BARREL PLATE WITH DCP HOLE 135* SHORT DIA BARREL 7 HOLE
288.108	GREENS DHS BARREL PLATE WITH DCP HOLE 135* SHORT DIA BARREL 8 HOLE

### Specification

- Plates available in Stainless steel.
- Plates available with 2,3,4,5,6,7 & 8 holes.
- It accepts 4.5mm Cortical screws & pediatricdhs screw.

### Indication for use

- These plates are used for simple and multiple pertrochanteric and intertrochanteric fractures, fractures of the femoral neck, and certain subtrochanteric fractures in Pediatric patients.