


GREENS - TUBULAR EXTERNAL FIXATORS FOR TIBIA, FEMUR & HUMERUS



CAT. NO.SS	Greens- Tubular External Fixators for Tibia, Femur & Humerus
361.030	Single Pin Clamp (Dx.)
361.035	Double Pin Clamp (Ord.)
361.040	Double Pin Clamp (Dx.)
361.050	Open Clamp (Dx.)
361.060	Connecting Clamp (Dx.)
361.070	Twin Adjustable Clamp Straight
361.080	Twin Adjustable Clamp Curved
361.090	Tube to Tube Clamp
361.100	Transverse Clamp
361.240	Clamp for Forearm 2.5mm x 4.0 mm
361.350	Clamp for Forearm 3.5mm x 4.0 mm
361.400	Connecting Clamp 4.0 mm x 4.0 mm
362.006	Spanner 06/ 07 mm
362.008	Box Spanner 8mm
362.011	Spanner 11mm
362.012	Socket Wrench in S.S. 11mm
362.015	Compression/ Distraction Device Large
362.020	Schanz Screw Introducer
362.050	Tripple Trocar with Handle 50mm
362.070	Tripple Trocarwith Handle 100mm
367.109	Tubular Rod (11mm) x 50mm
367.110	Tubular Rod (11mm) x 100mm
367.111	Tubular Rod (11mm) x 150mm
367.112	Tubular Rod (11mm) x 200mm
367.113	Tubular Rod (11mm) x 250mm
367.114	Tubular Rod (11mm) x 300mm
367.115	Tubular Rod (11mm) x 350mm
367.116	Tubular Rod (11mm) x 400mm
367.117	Tubular Rod (11mm) x 450mm
367.118	Tubular Rod (11mm) x 500mm
367.211	Tubular Rod (12mm) x 150mm Fessa
367.212	Tubular Rod (12mm) x 200mm Fessa
367.311	Tubular Rod (18mm) x 150mm Fessa
367.312	Tubular Rod (18mm) x 200mm Fessa
367.800	Tubular Rod (8mm) x 100mm
367.801	Tubular Rod (8mm) x 150mm
367.802	Tubular Rod (8mm) x 200mm
367.803	Tubular Rod (8mm) x 250mm



GREENS - TUBULAR EXTERNAL FIXATORS FOR TIBIA, FEMUR & HUMERUS



367,804	Tubular Rod (8mm) x 300mm
367,805	Tubular Rod (8mm) x 350mm
367,806	Tubular Rod (8mm) x 400mm
367,807	Tubular Rod (8mm) x 450mm
367,808	Tubular Rod (8mm) x 500mm
368,400	Connecting Rod Ø 4.0mm x 75 mm
368,401	Connecting Rod Ø 4.0mm x 100 mm
368,402	Connecting Rod Ø 4.0mm x 125 mm
368,403	Connecting Rod Ø 4.0mm x 150mm
368,404	Connecting Rod Ø 4.0mm x 200 mm
368,405	Connecting Rod Ø 4.0mm x 225 mm
368,406	Connecting Rod Ø 4.0mm x 250 mm
368,407	Connecting Rod Ø 4.0mm x 300 mm
368,408	Connecting Rod Ø 4.0mm x 175 mm
368,501	Connecting Rod Ø 5mm x 125 mm
368,502	Connecting Rod Ø 5mm x 150 mm
368,503	Connecting Rod Ø 5mm x 175 mm
368,504	Connecting Rod Ø 5mm x 200 mm
368,505	Connecting Rod Ø 5mm x 225 mm
368,506	Connecting Rod Ø 5mm x 250 mm
368,516	Connecting Rod Ø 5mm x 60 mm
368,517	Connecting Rod Ø 5mm x 70 mm
368,519	Connecting Rod Ø 5mm x 90 mm
368,521	Connecting Rod Ø 5mm x 110 mm
368,522	Connecting Rod Ø 5mm x 120 mm
368,528	Connecting Rod Ø 5mm x 180 mm
368,530	Connecting Rod Ø 5mm x 220 mm
369,100	Twin Rod Clamp Small (for upper limbs 2.5, 3.5mm Schanz Screws)
369,101	Twin Rod Clamp MEDIUM (for upper limbs 4.0, 4.0mm Schanz Screws)
369,200	Twin Pipe Clamp Large (for lower limbs 4.5, 5.0mm Schanz Screws)

Specification

- Greens - Tubular External Fixators for Tibia, Femur & Humerus available in Stainless steel.
- Tubular rod available in 50-500mm length with 8,11,12 & 18mm Dia.
- Connecting rod available in 75-30mm length with 4 & 5mm Dia.

Indication for use

- Greens - Tubular External fixator is used in these circumstances to stabilize the fracture adequately till the soft tissue condition is good enough for internal fixation.

