



CAT. NO.	METAL	HERBURT SCREW
230.012	SS	Cannulated Herbert Screw-, L- 12 mm
230.014	SS	Cannulated Herbert Screw-, L- 14 mm
230.016	SS	Cannulated Herbert Screw-, L- 16 mm
230.018	SS	Cannulated Herbert Screw-, L- 18 mm
230.020	SS	Cannulated Herbert Screw-, L- 20 mm
230.022	SS	Cannulated Herbert Screw-, L- 22 mm
230.024	SS	Cannulated Herbert Screw-, L- 24 mm
230.026	SS	Cannulated Herbert Screw-, L- 26 mm
230.028	SS	Cannulated Herbert Screw-, L- 28 mm
230.03	SS	Cannulated Herbert Screw-, L- 30 mm
230.212	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-12 MM
230.214	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-14 MM
230.216	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-16 MM
230.218	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-18 MM
230.22	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-20 MM
230.222	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-22 MM
230.224	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-24 MM
230.226	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-26 MM
230.228	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-28 MM
230.230	SS	2.4mm Headless compression Screw cannulated (HERBURT),L-30 MM
230.312	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-12 MM
230.314	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-14 MM
230.316	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-16 MM
230.318	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-18MM
230.320	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-20 MM
230.322	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-22 MM
230.324	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-24 MM
230.326	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-26 MM
230.328	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-28 MM
230.330	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-30 MM
230.332	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-32 MM

230.334	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-34 MM
230.336	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-36 MM
230.338	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-38 MM
230.340	SS	3.0mm Headless compression Screw cannulated(HERBURT),L-40 MM
639.012	TT	Cannulated Herbut Screw-, L- 12 mm
639.014	TT	Cannulated Herbut Screw-, L- 14 mm
639.016	TT	Cannulated Herbut Screw-, L- 16 mm
639.018	TT	Cannulated Herbut Screw-, L- 18 mm
639.02	TT	Cannulated Herbut Screw-, L- 20 mm
639.022	TT	Cannulated Herbut Screw-, L- 22 mm
639.024	TT	Cannulated Herbut Screw-, L- 24 mm
639.026	TT	Cannulated Herbut Screw-, L- 26 mm
639.028	TT	Cannulated Herbut Screw-, L- 28 mm
639.03	TT	Cannulated Herbut Screw-, L- 30 mm
639.032	TT	Cannulated Herbut Screw-, L- 32 mm
639.034	TT	Cannulated Herbut Screw-, L- 34 mm
639.035	TT	Cannulated Herbut Screw-, L- 35 mm
639.036	TT	Cannulated Herbut Screw-, L- 36 mm
639.038	TT	Cannulated Herbut Screw-, L- 38 mm
639.04	TT	Cannulated Herbut Screw-, L- 40 mm
639.042	TT	Cannulated Herbut Screw-, L- 42 mm
639.044	TT	Cannulated Herbut Screw-, L- 44 mm
639.045	TT	Cannulated Herbut Screw-, L- 45 mm
639.046	TT	Cannulated Herbut Screw-, L- 46 mm
639.048	TT	Cannulated Herbut Screw-, L- 48 mm
639.05	TT	Cannulated Herbut Screw-, L- 50 mm
639.055	TT	Cannulated Herbut Screw-, L- 55 mm
639.06	TT	Cannulated Herbut Screw-, L- 60 mm
639.065	TT	Cannulated Herbut Screw-, L- 65 mm
639.070	TT	Cannulated Herbut Screw-, L- 70 mm
639.210	TT	2.4mm Headless compression Screw cannulated (HERBURT),L-10 MM
639.212	TT	2.4mm Headless compression Screw cannulated (HERBURT),L-12 MM
639.214	TT	2.4mm Headless compression Screw cannulated (HERBURT),L-14 MM
639.216	TT	2.4mm Headless compression Screw cannulated (HERBURT),L-16 MM
639.218	TT	2.4mm Headless compression Screw cannulated (HERBURT),L-18 MM



639.220	TT	2.4mm Headless compression Screw cannulated (HERBURT),L -20 MM
639.222	TT	2.4mm Headless compression Screw cannulated (HERBURT),L -22 MM
639.224	TT	2.4mm Headless compression Screw cannulated (HERBURT),L -24 MM
639.226	TT	2.4mm Headless compression Screw cannulated (HERBURT),L -26 MM
639.228	TT	2.4mm Headless compression Screw cannula ted (HERBURT),L-28 MM
639.230	TT	2.4mm Headless compression Screw cannulated (HERBURT),L -30 MM
639.312	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -12 MM
639.314	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -14 MM
639.316	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -16 MM
639.318	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -18MM
639.320	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -20 MM
639.322	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -22 MM
639.324	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -24 MM
639.326	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -26 MM
639.328	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -28 MM
639.330	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -30 MM
639.332	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -32 MM
639.334	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -34 MM
639.336	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -36 MM
639.338	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -38 MM
639.340	TT	3.0mm Headless compression Screw cannulated(HERBURT),L -40 MM

Specification

- Available in 12-70mm length.
- Available in both Titanium and stainless steel.
- Available in Dia 3.5mm,2.4mm,3.0mm.

Indication for use

- A Herbur screw is a special purpose screw used to compress small fractures fragments. Herbert Screw Fixation of Scaphoid Frx. screw is threaded at both ends and can be countersunk beneath articular surface, obviating the need for removal after fix has healed.

