



CAT. NO.	METAL	High Tibial Osteotomies Locking Plate
643.833	TT	High Tibial Osteotomies Locking Plate 3H (Titanium)
643.834	TT	High Tibial Osteotomies Locking Plate 4H (Titanium)
643.835	TT	High Tibial Osteotomies Locking Plate 5H (Titanium)
643.836	TT	High Tibial Osteotomies Locking Plate 6H (Titanium)

Specification

- Plates available in Titanium material.
- Plate features locking compression holes and round locking holes that accepts 4.5/5.0mm cortical and Wise-lock screws.
- Plate available with 3,4,5 & 6 holes.

Indication for use

- It is indicated to fix Open-wedge and closed-wedge osteotomies of the Medial high tibia