

**HIP PLATES FOR INTERTROCHANTERIC
FEMORAL OSTEOTOMIES- 90°**



CAT. NO. SS	Hip Plates for Intertrochanteric Femoral Osteotomies- 90°
301.150	Hip Plates for Intertrochanteric Plates 90°, 4 Holes, Blade L- 50 mm
301.160	Hip Plates for Intertrochanteric Plates 90°, 4 Holes, Blade L- 60 mm
301.170	Hip Plates for Intertrochanteric Plates 90°, 4 Holes, Blade L- 70 mm

Specification

- Plates available in Stainless steel.
- Plates available in 4 holes with length 50,60 & 70mm.

Indication for use

- Intertrochanteric osteotomy alters the biomechanical environment of the fracture site and enhances fracture union.

