



CAT. NO. LVM	CAT.NO. TT	INFRAPECTINEAL PLATE
305.271L	700.991L	INFRAPECTINEAL PLATE -SMALL - LEFT
305.271R	700.991R	INFRAPECTINEAL PLATE -SMALL - RIGHT
305.275L	700.995L	INFRAPECTINEAL PLATE -LARGE - LEFT
305.275R	700.995R	INFRAPECTINEAL PLATE -LARGE - RIGHT

Specification

- Available in both Titanium & LVM.

Indication for use

- Infrapectineal plating provides rigid fixation for fracture patterns involving the quadrilateral surface; this is achieved by applying the plate in the same plane as the displacement rather than tangentially or perpendicular to the displacement plane.