



CAT. NO. SS	JEWETT NAIL PLATE 135*
287.402	JEWETT NAIL PLATE 135* 3 HOLE X 80MM
287.403	JEWETT NAIL PLATE 135* 3 HOLE X 85MM
287.404	JEWETT NAIL PLATE 135* 4 HOLE X 70MM
287.405	JEWETT NAIL PLATE 135* 4 HOLE X 75MM
287.406	JEWETT NAIL PLATE 135* 4 HOLE X 80MM
287.407	JEWETT NAIL PLATE 135* 4 HOLE X 85MM
287.414	JEWETT NAIL PLATE 135* 5 HOLE X 70MM
287.415	JEWETT NAIL PLATE 135* 5 HOLE X 75MM
287.416	JEWETT NAIL PLATE 135* 5 HOLE X 80MM
287.417	JEWETT NAIL PLATE 135* 5 HOLE X 85MM
287.421	JEWETT NAIL PLATE 135* 6 HOLE X 70MM
287.422	JEWETT NAIL PLATE 135* 6 HOLE X 75MM
287.423	JEWETT NAIL PLATE 135* 6 HOLE X 80MM
287.424	JEWETT NAIL PLATE 135* 6 HOLE X 85MM
287.426	JEWETT NAIL PLATE 135* 7 HOLE X 80MM
287.428	JEWETT NAIL PLATE 135* 8 HOLE X 75MM
287.430	JEWETT NAIL PLATE 135* 8 HOLE X 85MM
287.435	JEWETT NAIL PLATE 135* 9 HOLE X 80MM
287.456	JEWETT NAIL PLATE 135* 10 HOLE X 70MM
287.457	JEWETT NAIL PLATE 135* 10 HOLE X 75MM
287.458	JEWETT NAIL PLATE 135* 10 HOLE X 80MM
287.459	JEWETT NAIL PLATE 135* 10 HOLE X 85MM
287.463	JEWETT NAIL PLATE 135* 11 HOLE X 75MM
287.464	JEWETT NAIL PLATE 135* 11 HOLE X 80MM
287.470	JEWETT NAIL PLATE 135* 12 HOLE X 70MM
287.471	JEWETT NAIL PLATE 135* 12 HOLE X 75MM
287.472	JEWETT NAIL PLATE 135* 12 HOLE X 80MM
287.473	JEWETT NAIL PLATE 135* 12 HOLE X 85MM
287.477	JEWETT NAIL PLATE 135* 13 HOLE X 75MM
287.727	JEWETT NAIL PLATE 135* 8 HOLE X 70MM
287.729	JEWETT NAIL PLATE 135* 8 HOLE X 80MM

Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5,6,7,8,9,10,11,12 & 13 holes.
- Plate accepts 4.5mm Cortical screw.

Indication for use

- These plates are indicated for pertrochanteric, intertrochanteric, and basilar neck fractures.

