



CAT. NO. SS	Kirschner's Wire (K. Wire)(Sterilized)
241.103S	K. Wire Ø 1.0 mm - L 150mm(Sterilized)
241.121S	K. Wire Ø 1.2 mm - L 150mm(Sterilized)
241.151S	K. Wire Ø 1.5 mm - L 150mm(Sterilized)
241.152S	K. Wire Ø 1.5 mm - L 200mm(Sterilized)
241.181S	K. Wire Ø 1.8 mm - L 150mm(Sterilized)
241.201S	K. Wire Ø 2.0 mm - L 150mm(Sterilized)
241.202S	K. Wire Ø 2.0 mm - L 175mm(Sterilized)
241.203S	K. Wire Ø 2.0 mm - L 200mm(Sterilized)
241.251S	K. Wire Ø 2.5 mm - L 150mm(Sterilized)
241.252S	K. Wire Ø 2.5 mm - L 200mm(Sterilized)
241.301S	K. Wire Ø 3.0 mm - L 150mm(Sterilized)
241.302S	K. Wire Ø 3.0 mm - L 200mm(Sterilized)
241.421S	K. Wire Ø 1.5 mm - L 175mm(Sterilized)
241.426S	K. Wire Ø 2.5 mm - L 175mm(Sterilized)
241.431S	K. Wire Ø 3.0mm - L 175mm(Sterilized)

Specification

- Kirschner's Wire (K.wire) (Sterilised) available in Stainless steel & Titanium.
- Available length 150-175 with 1.0,1.2,1.5,1.8,2.0,2.5 & 3.0mm dia.

Indication for use

- Kirschner's Wire (K.wire) (Sterilised) are used for temporary fixation during some operations. After definitive fixation they are then removed. The pins are usually removed four to five weeks post operation