

L- BUTTRESS PLATES FOR 4.5/ 6.5 MM SCREWS



CAT. NO. SS	L- Buttress Plates for 4.5/ 6.5 mm Screws
298.103L	L- Buttress Plates - 03 Holes - Left, L- 69 mm
298.103R	L- Buttress Plates - 03 Holes - Right, L- 69 mm
298.104L	L- Buttress Plates - 04 Holes - Left, L- 85 mm
298.104R	L- Buttress Plates - 04 Holes - Right, L- 85 mm
298.105L	L- Buttress Plates - 05 Holes - Left, L- 101mm
298.105R	L- Buttress Plates - 05 Holes - Right, L- 101mm
298.106L	L- Buttress Plates - 06 Holes - Left, L- 117 mm
298.106R	L- Buttress Plates - 06 Holes - Right, L- 117 mm
298.107L	L- Buttress Plates - 07 Holes - Left, L- 133mm
298.107R	L- Buttress Plates - 07 Holes - Right, L- 133mm
298.108L	L- Buttress Plates - 08 Holes - Left, L- 149mm
298.108R	L- Buttress Plates - 08 Holes - Right, L- 149mm
298.110L	L- Buttress Plates - 10 Holes - Left,
298.110R	L- Buttress Plates - 10 Holes - Right,
298.112L	L- Buttress Plates - 12 Holes - Left,
298.112R	L- Buttress Plates - 12 Holes - Right,

Specification

- Plates available in Stainless steel.
- Plates available with 3,4,5,6,7,8,10 & 12 Holes.
- Plates accepts 4.5/6.5mm cortical screws and cancellous screw.

Indication for use

- L-Buttress Plate for the fracture of the tibia With double bend, as buttress for tibial and humeral head fractures.