



CAT. NO. SS	LC-DCP BROAD FOR 4.5 MM SCREWS
285.306	LC-DCP BROAD FOR4.5 MMSCREWS - 06 HOLES, L - 106 MM
285.307	LC-DCP BROAD FOR4.5 MMSCREWS - 07 HOLES, L - 124 MM
285.308	LC-DCP BROAD FOR4.5 MMSCREWS - 08 HOLES, L - 142 MM
285.309	LC-DCP BROAD FOR4.5 MMSCREWS - 09 HOLES, L - 160 MM
285.310	LC-DCP BROAD FOR4.5 MMSCREWS - 10 HOLES, L - 178 MM
285.311	LC-DCP BROAD FOR4.5 MMSCREWS - 11 HOLES, L - 196 MM
285.312	LC-DCP BROAD FOR4.5 MMSCREWS - 12 HOLES, L - 214 MM
285.313	LC-DCP BROAD FOR4.5 MMSCREWS - 13 HOLES,
285.314	LC-DCP BROAD FOR4.5 MMSCREWS - 14 HOLES, L - 250 MM
285.315	LC-DCP BROAD FOR4.5 MMSCREWS - 15 HOLES,
285.316	LC-DCP BROAD FOR4.5 MMSCREWS - 16 HOLES, L - 286 MM

Specification

- Available in Stainless Steel Material.
- Plate holes accept 4.5/6mm cortical screws and cancellous screws
- Plates available with 6,7,8,9,10,11,12,13,14,15 & 16 holes.
- LC-DCP holes allow compression in both the direction of the longitudinal axis
- Undercuts are designed for vascular preservation.

Indication for use

4.5mm Broad Dynamic Compression Plate indicated for fixation of-

- various long bones, such as the humerus, tibia and femur
- peri-prosthetic fractures and fractures in osteopenic bone
- nonunions or malunions in adult patients