



| CAT. NO.<br>SS | CAT. NO.<br>TT | LC-DCP SMALL FOR 3.5 MM SCREWS                       |
|----------------|----------------|--|
| 285.105        | 700.225        | LC-DCP SMALL FOR 3.5 MM SCREWS - 05 HOLES, L- 64 MM  |
| 285.106        | 700.226        | LC-DCP SMALL FOR 3.5 MM SCREWS - 06 HOLES, L- 77 MM  |
| 285.107        | 700.227        | LC-DCP SMALL FOR 3.5 MM SCREWS - 07 HOLES, L- 90 MM  |
| 285.108        | 700.228        | LC-DCP SMALL FOR 3.5 MM SCREWS - 08 HOLES, L- 103 MM |
| 285.109        | 700.229        | LC-DCP SMALL FOR 3.5 MM SCREWS - 09 HOLES, L- 116 MM |
| 285.110        | 700.230        | LC-DCP SMALL FOR 3.5 MM SCREWS - 10 HOLES, L- 129 MM |
| 285.111        | 700.231        | LC-DCP SMALL FOR 3.5 MM SCREWS - 11 HOLES, L- 142 MM |
| 285.112        | 700.232        | LC-DCP SMALL FOR 3.5 MM SCREWS - 12 HOLES, L- 155 MM |

## Specification

- Available in both Stainless Steel and Titanium material.
- The oval shape holes are designed to create dynamic pressure between the fractured bone fragments.
- All holes of the plate support 3.5mm cortical & 4.0mm cancellous screws.

## Indication for use

- It is used to fix fractures and non-unions of the fibula, Distal tibia, Pelvis, Ulna, Radius, Humerous, Olecranon, Scapula & Clavicle.

