



CAT. NO. SS	CAT. NO. TT	LC-LCP CLAVICLE HOOK PLATE
321.116L	651.104L	LC-LCP CLAVICLE HOOK PLATE - LEFT -4H
321.116R	651.104R	LC-LCP CLAVICLE HOOK PLATE - RIGHT -4H
321.117L	651.105L	LC-LCP CLAVICLE HOOK PLATE - LEFT -5H
321.117R	651.105R	LC-LCP CLAVICLE HOOK PLATE - RIGHT -5H
321.118L	651.106L	LC-LCP CLAVICLE HOOK PLATE - LEFT -6H
321.118R	651.106R	LC-LCP CLAVICLE HOOK PLATE - RIGHT -6H
321.119L	651.107L	LC-LCP CLAVICLE HOOK PLATE - LEFT -7H
321.119R	651.107R	LC-LCP CLAVICLE HOOK PLATE - RIGHT -7H

Specification

- Dynamic compression screw holes accept 3.5mm Cortical & locking Screws and 4.0mm cancellous bone screws.
- Plate is compatible with the 3.5mm DCP drill guide.
- Anteriolateral screw hole provides additional option for screw fixation in the lateral clavicle.
- Hooks provide additional support to both lateral clavicle fractures and acromioclavicular joint dislocations.
- Plates available with 4, 5, 6 or 7 holes.
- Offset hook design to avoid insertion of hook into acromioclavicular ligament.

Indication for use

The 3.5mm clavicle hook plate is indicated for -

- Lateral clavicle fracture
- Dislocations of the acromioclavicular joint