



CAT. NO. SS	CAT. NO. TT	LC-LCP SMALL FOR Ø 3.5MM
310.005	610.005	LC-LCP SMALL FOR Ø 3.5MM- 5H
310.006	610.006	LC-LCP SMALL FOR Ø 3.5MM- 6H
310.007	610.007	LC-LCP SMALL FOR Ø 3.5MM- 7H
310.008	610.008	LC-LCP SMALL FOR Ø 3.5MM- 8H
310.009	610.009	LC-LCP SMALL FOR Ø 3.5MM- 9H
310.010	610.010	LC-LCP SMALL FOR Ø 3.5MM- 10H
310.011	610.011	LC-LCP SMALL FOR Ø 3.5MM- 11H
310.012	610.012	LC-LCP SMALL FOR Ø 3.5MM- 12H
310.013	610.013	LC-LCP SMALL FOR Ø 3.5MM- 13H
310.014	610.014	LC-LCP SMALL FOR Ø 3.5MM- 14H

**Specification**

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plates available with 5,6,7,8,9,10,11,12,13 & 14 holes.

**Indication for use**

- Intraarticular fractures of the distal humerus
- Supracondylar fractures of the distal humerus
- Non-unions of the distal humerus