



| CAT. NO. SS | CAT. NO. TT | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM |
|----------------|----------------|--|
| 313.103 | 613.103 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 03 HOLES |
| 313.104 | 613.104 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 04 HOLES |
| 313.105 | 613.105 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 05 HOLES |
| 313.106 | 613.106 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 06 HOLES |
| 313.107 | 613.107 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 07 HOLES |
| 313.108 | 613.108 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 08 HOLES |
| 313.109 | 613.109 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 09 HOLES |
| 313.110 | 613.110 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 10 HOLES |
| 313.111 | 613.111 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 11 HOLES |
| 313.112 | 613.112 | LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 12 HOLES |

Specification

- One-Third Tubular plate is available in both Stainless Steel and Titanium material.
- The plate holes support 3.5mm cortex and 4.0mm cancellous screws.

Indication for use

- It is used as buttress/supporting plate for tubular Metatarsals and Metacarpals.