



CAT. NO. SS	CAT. NO. TT	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM
313.103	613.103	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 03 HOLES
313.104	613.104	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 04 HOLES
313.105	613.105	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 05 HOLES
313.106	613.106	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 06 HOLES
313.107	613.107	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 07 HOLES
313.108	613.108	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 08 HOLES
313.109	613.109	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 09 HOLES
313.110	613.110	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 10 HOLES
313.111	613.111	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 11 HOLES
313.112	613.112	LC-LCP 1/3RD TUBULAR PLATES FOR 3.5 MM- 12 HOLES

Specification

- One-Third Tubular plate is available in both Stainless Steel and Titanium material.
- The plate holes support 3.5mm cortex and 4.0mm cancellous screws.

Indication for use

- It is used as buttress/supporting plate for tubular Metatarsals and Metacarpals.

