



CAT. NO. SS	CAT. NO. TT	LC-LCP BROAD FOR Ø 4.5MM
310.206	610.206	LC-LCP BROAD FOR Ø 4.5MM - 6H
310.207	610.207	LC-LCP BROAD FOR Ø 4.5MM - 7H
310.208	610.208	LC-LCP BROAD FOR Ø 4.5MM - 8H
310.209	610.209	LC-LCP BROAD FOR Ø 4.5MM - 9H
310.210	610.210	LC-LCP BROAD FOR Ø 4.5MM - 10H
310.211	610.211	LC-LCP BROAD FOR Ø 4.5MM - 11H
310.212	610.212	LC-LCP BROAD FOR Ø 4.5MM - 12H
310.213	610.213	LC-LCP BROAD FOR Ø 4.5MM - 13H
310.214	610.214	LC-LCP BROAD FOR Ø 4.5MM - 14H
310.215	610.215	LC-LCP BROAD FOR Ø 4.5MM - 15H
310.216	610.216	LC-LCP BROAD FOR Ø 4.5MM - 16H
310.217	610.217	LC-LCP BROAD FOR Ø 4.5MM - 17H
310.218	610.218	LC-LCP BROAD FOR Ø 4.5MM - 18H

Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm/5.0mm cortical and locking screws.
- Plate available with 6,7,8,9,10,11,12,13,14,15,16,17 & 18 holes.

Indication for use

- It is indicated for fixation of various long bones, such as the humerus, femur and tibia and for use in fixation of peri-prosthetic fractures, osteopenic bone and fixation of non-unions or malunions in adult patients.