



CAT. NO. SS	CAT. NO. TT	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM
333.030	666.530	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 30MM
333.035	666.535	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 35MM
333.040	666.540	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 40MM
333.045	666.545	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 45MM
333.050	666.550	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 50MM
333.055	666.555	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 55MM
333.060	666.560	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 60MM
333.065	666.565	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 65MM
333.070	666.570	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 70MM
333.075	666.575	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 75MM
333.080	666.580	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 80MM
333.085	666.585	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 85MM
333.090	666.590	LC-LCP CANCELLOUS (16MM) SCREW Ø 5.0MM X 90MM

Specification

- Available in 30mm to 90mm length with 5.0mm diameter
- Partially 16mm threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones