



CAT. NO. SS	CAT. NO. TT	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM
334.030	634.030	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 30MM
334.035	634.035	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 35MM
334.040	634.040	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 40MM
334.045	634.045	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 45MM
334.050	634.050	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 50MM
334.055	634.055	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 55MM
334.060	634.060	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 60MM
334.065	634.065	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 65MM
334.070	634.070	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 70MM
334.075	634.075	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 75MM
334.080	634.080	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 80MM
334.085	634.085	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 85MM
334.090	634.090	LC-LCP CANCELLOUS (32MM) SCREW Ø 5.0MM X 90MM

Specification

- Available in 30mm to 90mm length with 5.0mm diameter
- Partially 32mm threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhances the holding power of the screw in soft bones