



CAT. NO. SS	LC-LCP Cloverleaf Plate for Ø 3.5mm
317.305	LC-LCP Cloverleaf Plate for Ø 3.5mm - 5H
317.30 6	LC-LCP Cloverleaf Plate for Ø 3.5mm - 6H
317.30 7	LC-LCP Cloverleaf Plate for Ø 3.5mm - 7H
317.30 8	LC-LCP Cloverleaf Plate for Ø 3.5mm - 8H

## Specification

- Plates available in Stainless steel.
- Plate accepts 3.5mm cortical and Lcp 3.5mm screws.
- Plate available with 5,6,7 & 8 holes.

## Indication for use

- These plates are indicated for the fractures of Proximal Humerus and Distal Tibia. One or both shoulders of the Cloverleaf Plate can be cut off with large Wire Cutter.