



CAT. NO. SS	CAT. NO. TT	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM
316.403L	616.403L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 3H- LEFT
316.403R	616.403R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 3H- RIGHT
316.404L	616.404L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 4H- LEFT
316.404R	616.404R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 4H- RIGHT
316.405L	616.405L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 5H- LEFT
316.405R	616.405R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 5H- RIGHT
316.406L	616.406L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 6H- LEFT
316.406R	616.406R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 6H- RIGHT
316.407L	616.407L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 7H- LEFT
316.407R	616.407R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 7H- RIGHT
316.408L	616.408L	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 8H- LEFT
316.408R	616.408R	LC-LCP L.BUTTRESS PLATE FOR Ø 4.5MM - 8H- RIGHT

## Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5/5.0mm cortical and Wise-lock screws.
- Plate available with 3,4,5,6,7 & 8 holes.

## Indication for use

- 4.5/5.0mm Wise Lock 'L' Buttress Plate is indicated for fixation of fractures, osteotomies and non-unions of the clavicle, scapula, olecranon, humerus, radius ulna, pelvis, distal tibia, fibula, particularly in osteopenic bones.