

## LC-LCP MEDIAL DISTAL HUMERAL PLATES 2.7/3.5MM



| CAT. NO.<br>SS | CAT. NO.<br>TT | LC-LCP MEDIAL DISTAL HUMERAL PLATES 2.7/3.5MM               |
|----------------|----------------|-------------------------------------------------------------|
| 448.503L       | 648.503L       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 3H- LEFT  |
| 448.503R       | 648.503R       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 3H- RIGHT |
| 448.505L       | 648.505L       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 5H- LEFT  |
| 448.505R       | 648.505R       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 5H- RIGHT |
| 448.507L       | 648.507L       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 7H- LEFT  |
| 448.507R       | 648.507R       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 7H- RIGHT |
| 448.509L       | 648.509L       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 9H- LEFT  |
| 448.509R       | 648.509R       | LC-LCP MEDIAL DISTAL HUMERAL PLATES Ø 2.7/3.5MM - 9H- RIGHT |

### Specification

- Plates available in Stainless steel & Titanium material.
- Plate accepts 2.7/3.5mm cortical and locking screws.
- Plate available with 3,5,7 & 9 Holes.
- Plates are left & right configured.

### Indication for use

- It is indicated for intra-articular fracture of the distal humerus, comminuted  
Supracondylar fractures, osteotomies and non-unions of the distal humerus.

