

LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM (WITHOUT TAB)



CAT. NO. SS	CAT. NO. TT	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM (WITHOUT TAB)
322.804L	622.804L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 4H- LEFT
322.804R	622.804R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 4H- RIGHT
322.806L	622.806L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 6H- LEFT
322.806R	622.806R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 6H- RIGHT
322.808L	622.808L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 8H- LEFT
322.808R	622.808R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM -(WITHOUT TAB) 8H - RIGHT
322.810L	622.810L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 10H- LEFT
322.810R	622.810R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 10H- RIGHT
322.812L	622.812L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 12H- LEFT
322.812R	622.812R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 12H- RIGHT
322.814L	622.814L	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 14H- LEFT
322.814R	622.814R	LC-LCP MEDIAL DISTAL TIBIA PLATE 3.5MM(WITHOUT TAB) - 14H- RIGHT

Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plate available with 4,6,8,10,12 & 14 holes

Indication for use

Wise Lock Medial Distal Tibia Plates without Tab are intended for :

- Fixation of complex intra- and extra-articular fractures of the distal tibia
- Osteotomies of the distal tibia