



| CAT. NO. SS | CAT. NO. TT | LC-LCP NARROW FOR Ø 4.5MM |
|----------------|----------------|---------------------------------|
| 310.106 | 610.106 | LC-LCP NARROW FOR Ø 4.5MM - 6H |
| 310.107 | 610.107 | LC-LCP NARROW FOR Ø 4.5MM - 7H |
| 310.108 | 610.108 | LC-LCP NARROW FOR Ø 4.5MM - 8H |
| 310.109 | 610.109 | LC-LCP NARROW FOR Ø 4.5MM - 9H |
| 310.110 | 610.110 | LC-LCP NARROW FOR Ø 4.5MM - 10H |
| 310.111 | 610.111 | LC-LCP NARROW FOR Ø 4.5MM - 11H |
| 310.112 | 610.112 | LC-LCP NARROW FOR Ø 4.5MM - 12H |
| 310.114 | 610.114 | LC-LCP NARROW FOR Ø 4.5MM - 14H |
| 310.116 | 610.116 | LC-LCP NARROW FOR Ø 4.5MM - 16H |

Specification

- Available in both Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5mm cortical and 5.0mm locking screws.
- Plate available with 6,7,8,9,10,11,12,14 &16 holes.

Indication for use

- Lc- Lcp Narrow Plate with LC under cuts 4.5mm is indicated for fixation of various long bones, such as the humerus, femur and tibia and for use in fixation of peri-prosthetic fractures, osteopenic bone and fixation of non-unions or malunions in adult patients.