



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 2.4 MM (HEXAGONAL HEAD)
325.010	627.010	LC - LCP SCREW 2.4 MM 10 MM (HEXAGONAL HEAD)
325.012	627.012	LC - LCP SCREW 2.4 MM 12 MM (HEXAGONAL HEAD)
325.014	627.014	LC - LCP SCREW 2.4 MM 14 MM (HEXAGONAL HEAD)
325.016	627.016	LC - LCP SCREW 2.4 MM 16 MM (HEXAGONAL HEAD)
325.018	627.018	LC - LCP SCREW 2.4 MM 18 MM (HEXAGONAL HEAD)
325.020	627.020	LC - LCP SCREW 2.4 MM 20 MM (HEXAGONAL HEAD)
325.022	627.022	LC - LCP SCREW 2.4 MM 22 MM (HEXAGONAL HEAD)
325.024	627.024	LC - LCP SCREW 2.4 MM 24 MM (HEXAGONAL HEAD)
325.026	627.026	LC - LCP SCREW 2.4 MM 26 MM (HEXAGONAL HEAD)
325.028	627.028	LC - LCP SCREW 2.4 MM 28 MM (HEXAGONAL HEAD)
325.030	627.030	LC - LCP SCREW 2.4 MM 30 MM (HEXAGONAL HEAD)
325.032	627.034	LC - LCP SCREW 2.4 MM 32 MM (HEXAGONAL HEAD)
325.034	627.032	LC - LCP SCREW 2.4 MM 34 MM (HEXAGONAL HEAD)
325.036	627.036	LC - LCP SCREW 2.4 MM 36 MM (HEXAGONAL HEAD)
325.038	627.038	LC - LCP SCREW 2.4 MM 38 MM (HEXAGONAL HEAD)
325.040	627.040	LC - LCP SCREW 2.4 MM 40 MM (HEXAGONAL HEAD)

Specification

- Compress the plate to the bone or create axial compression.
- Available in 10 mm to 40 mm length with 2.4 mm diameter.
- Fully threaded shaft.
- Available in Stainless steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Wise-Lock Plates.