



CAT. NO. SS	CAT. NO. TT	LCP VARIABLE ANGLE SCREW 2.4MM (STAR HEAD)
325.050	627.210	LCP VARIABLE ANGLE SCREW 2.4MM 10.0MM (STAR HEAD)
325.052	627.212	LCP VARIABLE ANGLE SCREW 2.4MM 12.0MM (STAR HEAD)
325.054	627.214	LCP VARIABLE ANGLE SCREW 2.4MM 14.0MM (STAR HEAD)
325.056	627.216	LCP VARIABLE ANGLE SCREW 2.4MM 16.0MM (STAR HEAD)
325.058	627.218	LCP VARIABLE ANGLE SCREW 2.4MM 18.0MM (STAR HEAD)
325.060	627.220	LCP VARIABLE ANGLE SCREW 2.4MM 20.0MM (STAR HEAD)
325.062	627.222	LCP VARIABLE ANGLE SCREW 2.4MM 22.0MM (STAR HEAD)
325.064	627.224	LCP VARIABLE ANGLE SCREW 2.4MM 24.0MM (STAR HEAD)
325.066	627.226	LCP VARIABLE ANGLE SCREW 2.4MM 26.0MM (STAR HEAD)
325.068	627.228	LCP VARIABLE ANGLE SCREW 2.4MM 28.0MM (STAR HEAD)
325.070	627.230	LCP VARIABLE ANGLE SCREW 2.4MM 30.0MM (STAR HEAD)
325.072	627.232	LCP VARIABLE ANGLE SCREW 2.4MM 32.0MM (STAR HEAD)
325.074	627.234	LCP VARIABLE ANGLE SCREW 2.4MM 34.0MM (STAR HEAD)
325.076	627.236	LCP VARIABLE ANGLE SCREW 2.4MM 36.0MM (STAR HEAD)
325.078	627.238	LCP VARIABLE ANGLE SCREW 2.4MM 38.0MM (STAR HEAD)
325.080	627.240	LCP VARIABLE ANGLE SCREW 2.4MM 40.0MM (STAR HEAD)

Specification

- Self- tapping screw can tap its own hole and hence can ease the process of insertion.
- Available in 10mm to 40mm length with 2.4mm diameter.
- Fully threaded shaft.
- Available in both Stainless Steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Plates.

