



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 2.7 MM LENGTH STAR HEAD
326.310	626.310	LC - LCP SCREW 2.7 MM LENGTH 10 MM STAR HEAD
326.312	626.312	LC - LCP SCREW 2.7 MM LENGTH 12 MM STAR HEAD
326.314	626.314	LC - LCP SCREW 2.7 MM LENGTH 14 MM STAR HEAD
326.316	626.316	LC - LCP SCREW 2.7 MM LENGTH 16 MM STAR HEAD
326.318	626.318	LC - LCP SCREW 2.7 MM LENGTH 18 MM STAR HEAD
326.320	626.320	LC - LCP SCREW 2.7 MM LENGTH 20 MM STAR HEAD
326.322	626.322	LC - LCP SCREW 2.7 MM LENGTH 22 MM STAR HEAD
326.324	626.324	LC - LCP SCREW 2.7 MM LENGTH 24 MM STAR HEAD
326.326	626.326	LC - LCP SCREW 2.7 MM LENGTH 26 MM STAR HEAD
326.328	626.328	LC - LCP SCREW 2.7 MM LENGTH 28 MM STAR HEAD
326.330	626.330	LC - LCP SCREW 2.7 MM LENGTH 30 MM STAR HEAD
326.332	626.332	LC - LCP SCREW 2.7 MM LENGTH 32 MM STAR HEAD
326.334	626.334	LC - LCP SCREW 2.7 MM LENGTH 34 MM STAR HEAD
326.336	626.336	LC - LCP SCREW 2.7 MM LENGTH 36 MM STAR HEAD
326.338	626.338	LC - LCP SCREW 2.7 MM LENGTH 38 MM STAR HEAD
326.340	626.340	LC - LCP SCREW 2.7 MM LENGTH 40 MM STAR HEAD

Specification

- Self-tapping screw can tap its own hole and hence can ease the process of insertion.
- Available in 10mm to 40mm length with 2.7mm diameter.
- Fully threaded shaft.
- Available in both Stainless Steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Plates.