

LC - LCP SCREW 3.5 MM LENGTH - STAR HEAD



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 3.5 MM LENGTH - STAR HEAD
326.512	626.512	LC - LCP SCREW 3.5 MM LENGTH 12 MM - STAR HEAD
326.514	626.514	LC - LCP SCREW 3.5 MM LENGTH 14 MM - STAR HEAD
326.516	626.516	LC - LCP SCREW 3.5 MM LENGTH 16 MM - STAR HEAD
326.518	626.518	LC - LCP SCREW 3.5 MM LENGTH 18 MM - STAR HEAD
326.520	626.520	LC - LCP SCREW 3.5 MM LENGTH 20 MM - STAR HEAD
326.522	626.522	LC - LCP SCREW 3.5 MM LENGTH 22 MM - STAR HEAD
326.524	626.524	LC - LCP SCREW 3.5 MM LENGTH 24 MM - STAR HEAD
326.526	626.526	LC - LCP SCREW 3.5 MM LENGTH 26 MM - STAR HEAD
326.528	626.528	LC - LCP SCREW 3.5 MM LENGTH 28 MM - STAR HEAD
326.530	626.530	LC - LCP SCREW 3.5 MM LENGTH 30 MM - STAR HEAD
326.535	626.535	LC - LCP SCREW 3.5 MM LENGTH 35 MM - STAR HEAD
326.540	626.540	LC - LCP SCREW 3.5 MM LENGTH 40 MM - STAR HEAD
326.545	626.545	LC - LCP SCREW 3.5 MM LENGTH 45 MM - STAR HEAD
326.550	626.550	LC - LCP SCREW 3.5 MM LENGTH 50 MM - STAR HEAD
326.555	626.555	LC - LCP SCREW 3.5 MM LENGTH 55 MM - STAR HEAD
326.560	626.560	LC - LCP SCREW 3.5 MM LENGTH 60 MM - STAR HEAD
326.565	626.565	LC - LCP SCREW 3.5 MM LENGTH 65 MM - STAR HEAD
326.570	626.570	LC - LCP SCREW 3.5 MM LENGTH 70 MM - STAR HEAD
326.575	626.575	LC - LCP SCREW 3.5 MM LENGTH 75 MM - STAR HEAD
326.580	626.580	LC - LCP SCREW 3.5 MM LENGTH 80 MM - STAR HEAD

Specification

- Self-tapping screw can tap its own hole and hence can ease the process of insertion.
- Available in 12mm to 80mm length with 3.5mm diameter.
- Fully threaded shaft.
- Available in both Stainless Steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Plates.