



CAT. NO. SS	CAT. NO. TT	LC - LCP SCREW 5.0 MM LENGTH - STAR HEAD
326.612	626.612	LC - LCP SCREW 5.0 MM LENGTH 12 MM - STAR HEAD
326.614	626.614	LC - LCP SCREW 5.0 MM LENGTH 14 MM - STAR HEAD
326.616	626.616	LC - LCP SCREW 5.0 MM LENGTH 16 MM - STAR HEAD
326.618	626.618	LC - LCP SCREW 5.0 MM LENGTH 18 MM - STAR HEAD
326.620	626.620	LC - LCP SCREW 5.0 MM LENGTH 20 MM - STAR HEAD
326.622	626.622	LC - LCP SCREW 5.0 MM LENGTH 22 MM - STAR HEAD
326.624	626.624	LC - LCP SCREW 5.0 MM LENGTH 24 MM - STAR HEAD
326.626	626.626	LC - LCP SCREW 5.0 MM LENGTH 26 MM - STAR HEAD
326.628	626.628	LC - LCP SCREW 5.0 MM LENGTH 28 MM - STAR HEAD
326.630	626.630	LC - LCP SCREW 5.0 MM LENGTH 30 MM - STAR HEAD
326.635	626.635	LC - LCP SCREW 5.0 MM LENGTH 35 MM - STAR HEAD
326.640	626.640	LC - LCP SCREW 5.0 MM LENGTH 40 MM - STAR HEAD
326.645	626.645	LC - LCP SCREW 5.0 MM LENGTH 45 MM - STAR HEAD
326.650	626.650	LC - LCP SCREW 5.0 MM LENGTH 50 MM - STAR HEAD
326.655	626.655	LC - LCP SCREW 5.0 MM LENGTH 55 MM - STAR HEAD
326.660	626.660	LC - LCP SCREW 5.0 MM LENGTH 60 MM - STAR HEAD
326.665	626.665	LC - LCP SCREW 5.0 MM LENGTH 65 MM - STAR HEAD
326.670	626.670	LC - LCP SCREW 5.0 MM LENGTH 70 MM - STAR HEAD
326.675	626.675	LC - LCP SCREW 5.0 MM LENGTH 75 MM - STAR HEAD
326.680	626.680	LC - LCP SCREW 5.0 MM LENGTH 80 MM - STAR HEAD

Specification

- Self- tapping screw can tap its own hole and hence can ease the process of insertion.
- Available in 12mm to 80mm length with 5.0mm diameter.
- Fully threaded shaft.
- Available in both Stainless Steel & Titanium.

Indication for use

- It is used in the fixation of Small Fragment Plates.