



| CAT. NO. SS | CAT. NO. TT | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM |
|----------------|----------------|--|
| 315.404 | 615.404 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 4H |
| 315.405 | 615.405 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 5H |
| 315.406 | 615.406 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 6H |
| 315.407 | 615.407 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 7H |
| 315.408 | 615.408 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 8H |
| 315.409 | 615.409 | LC-LCP T.BUTTRESS PLATE FOR Ø 4.5MM - 9H |

Specification

- Available in both Stainless Steel and Titanium material
- Plate holes accept 4.5mm cortical & locking 5.0mm screws
- Plate available with 4,5,6,7,8 & 9 holes

Indication for use

- Lateral plateau fracture of the tibia.
- Areas of thin metaphyseal cortical bone or to prevent the collapsing during the treatment of cancellous defects in that area.

