



CAT. NO. SS	CAT. NO. TT	LC-LCP T.Plate for Ø 4.5 / 5.0mm
314.403	614.403	LC-LCP TPLATE FOR Ø 4.5MM- 3H
314.404	614.404	LC-LCP TPLATE FOR Ø 4.5MM- 4H
314.405	614.405	LC-LCP TPLATE FOR Ø 4.5MM- 5H
314.406	614.406	LC-LCP TPLATE FOR Ø 4.5MM- 6H
314.408	614.408	LC-LCP TPLATE FOR Ø 4.5MM- 8H
314.409	614.409	LC-LCP TPLATE FOR Ø 4.5MM- 9H
314.410	614.410	LC-LCP TPLATE FOR Ø 4.5MM- 10H

Specification

- Plates (Right & Left) available in both-Titanium & Stainless steel.
- Plate accepts 4.5/5.0mm cortical and Wise-lock screws.
- Plates are available with 3,4,5,6,8,9 & 10 holes.

Indication for use

- Lc-Lcp "T" Plate 3.5mm is used for fixation of Partial Articular or Split Metaphyseal fractures of long bones, Proximal and Distal part of Humerus, Radius, Ulna, Tibia and Fibula.