




CAT.NO. SS	CAT. NO. TT	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM
320.405	620.405	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 5H- LEFT
320.407	620.407	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 7H- LEFT
320.409	620.409	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 9H- LEFT
320.411	620.411	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 11H- LEFT
320.413	620.413	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 13H- LEFT
320.415	620.415	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 15H- LEFT
320.505	620.505	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 5H- RIGHT
320.507	620.507	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 7H- RIGHT
320.509	620.509	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 9H- RIGHT
320.511	620.511	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 11H- RIGHT
320.513	620.513	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 13H- RIGHT
320.515	620.515	LC-LCP DISTAL FEMORAL PLATE FOR Ø 4.5MM - 15H- RIGHT

Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5/5.0mm cortical and locking screws.
- Plate available with 5,7,9,11,13 & 15 holes.

Indication for use

- Distal shaft fractures
- Supracondylar fractures
- Intra-articular fractures
- Periprosthetic fractures