



CAT. NO. TT	LCP Condylar Plate 2.0 MM
604.167	LCP Condylar Plate 2.0 - 4H
604.170	LCP Condylar Plate 2.0 - 7H

## Specification

- Plates available in Titanium.
- Condylar Plate allow implant placement to address the individual fracture pattern.
- Plates available with 4 & 7 holes.

## Indication for use

- Locking Adaption Plate is intended for temporary fixation, correction or stabilization of bones in anatomical region of the hand.

### It is used in

- Fracture fixation of the phalanges and metacarpals
- Osteotomies
- Arthrodeses
- Replantations and reconstructions of phalanges and metacarpals, particularly in osteopenic bone

