



CAT. NO. SS	CAT. NO. TT	LCP DISTAL RADIUS STRAIGHT PLATE	
350.273	644.273	LCP DISTAL RADIUS STRAIGHT PLATE	- DORSAL PLATES 5HOLES
350.274	644.274	LCP DISTAL RADIUS STRAIGHT PLATE	- DORSAL PLATES 6HOLES

## Specification

- Plates available in Stainless steel & Titanium material.
- Plate features locking compression holes and round locking holes that accepts 2.4mm cortical and locking screws.
- Plate available with 5 & 6 holes
- Locking Plates provides anatomical reduction, i.e., adaptation of fracture fragments to restore normal anatomy.
- It gives stability in fixation to the fracture. Plate provides optimal compression and stability.
- Plate preserves blood supply to the soft tissues and bone.

## Indication for use

- For the fixation of complete intra- and extra-articular fractures and osteotomies of the distal radius extending to diaphysis and metaphysis.