



CAT. NO. SS	CAT. NO. TT	LCP DISTAL RADIUS T - PLATE - DORSAL PLATES
349.273	649.273	LCP DISTAL RADIUS T - PLATE - DORSAL PLATES, 3HOLES SHAFT,3HOLES HEAD
349.274	649.274	LCP DISTAL RADIUS T - PLATE - DORSAL PLATES, 4HOLES SHAFT,3HOLES HEAD

Specification

- Plates available in Stainless steel & Titanium material.
- Plate features locking compression holes and round locking holes that accepts 2.4mm cortical and Wise-lock screws.
- Plate available with 3 head holes & 3,4 shaft holes Locking Plates provides anatomical reduction, i.e., adaptation of fracture fragments to restore normal anatomy.
- It gives stability in fixation to the fracture. Plate provides optimal compression and stability.
- Plate preserves blood supply to the soft tissues and bone.
- Plate preserves blood supply to soft tissues and bone.

Indication for use

- These plates are indicated for extra-articular, partial articular & complete articular fractures of distal radius.