

**LCP SUPERIOR ANTERIOR CLAVICLE
PLATE 3.5MM**



Cat No. SS	Cat No. TT	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM
351.206L	651.206L	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 6 HOLE, LEFT
351.206R	651.206R	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 6 HOLE, RIGHT
351.207L	651.207L	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 7 HOLE, LEFT
351.207R	651.207R	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 7 HOLE, RIGHT
351.208L	651.208L	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 8 HOLE, LEFT
351.208R	651.208R	LCP SUPERIOR ANTERIOR CLAVICLE PLATE 3.5MM, 8 HOLE, RIGHT

Specification

- Pre-contoured plate for anatomical shape
- Available with 6, 7 and 8 holes (Left and Right)
- Plate features locking compression holes and round locking holes that accept 3.5mm cortical and locking screws
- Available in Stainless steel and Titanium

Indication for use

- LCP Clavicle Plate is used for the fixation of fractures, malunions, nonunions and osteotomies of the clavicle.