

# LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM



| Cat. No.<br>SS | Cat. No.<br>TT | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM                |
|----------------|----------------|----------------------------------------------------------------------------------|
| 352.103L       | 652.103L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 3 HOLE, LEFT  |
| 352.103R       | 652.103R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 3 HOLE, RIGHT |
| 352.104L       | 652.104L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 4 HOLE, LEFT  |
| 352.104R       | 652.104R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 4 HOLE, RIGHT |
| 352.105L       | 652.105L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 5 HOLE, LEFT  |
| 352.105R       | 652.105R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 5 HOLE, RIGHT |
| 352.106L       | 652.106L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 6 HOLE, LEFT  |
| 352.106R       | 652.106R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 6 HOLE, RIGHT |
| 352.107L       | 652.107L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 7 HOLE, LEFT  |
| 352.107R       | 652.107R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 7 HOLE, RIGHT |
| 352.108L       | 652.108L       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 8 HOLE, LEFT  |
| 352.108R       | 652.108R       | LCP SUPERIOR ANTERIOR CLAVICLE PLATE WITH LATERAL EXTENSION 3.5MM, 8 HOLE, RIGHT |

## Specification

- Pre-contoured plate for anatomical shape
- plates with lateral extension are available in six lengths with 3, 4, 5, 6, 7 or 8 Combi holes in the shaft
- Plates with lateral extension feature six distal locking holes that accept 2.7 mm locking and cortical screws. The Combi holes in the plate shaft accept 3.5 mm locking, 3.5 mm cortex, and 4.0 mm cancellous bone screws.
- Available in Stainless steel and Titanium, Left & Right

## Indication for use

- Plates fit on the anterior aspect of the clavicle
- Notches in the plate allow any necessary plate contouring
- The Combi holes allow fixation with locking screws in the threaded section for angular stability, and cortex screws in the dynamic compression unit section for compression.