



CAT. NO. SS	CAT. NO. TT	LCP CALCANEAL PLATE - PLAIN
318.060L	618.00 6L	LCP CALCANEAL PLATE, L - 60MM - LEFT
318.060R	618.006R	LCP CALCANEAL PLATE, L - 60MM -RIGHT
318.070L	618.00 7L	LCP CALCANEAL PLATE, L - 70MM - LEFT
318.070R	618.00 7R	LCP CALCANEAL PLATE, L - 70MM -RIGHT

Specification

- This plate are available in 60-70mm length.
- All holes of the plate support Lcp 3.5mm screws.
- Available in both Titanium & Stainless Steel.

Indication for use

- 3.5mm Calcaneal Plate is used for non-articular body fractures (Non-comminuted and comminuted).

