



CAT. NO. SS	CAT. NO. TT	LCP CALCANEUS PLATE - AO TYPE
318.160L	618.060L	LCP CALCANEUS PLATE, L - 60MM - LEFT
318.160R	618.060R	LCP CALCANEUS PLATE, L - 60MM - RIGHT
318.170L	618.070L	LCP CALCANEUS PLATE, L - 70MM - LEFT
318.170R	618.070R	LCP CALCANEUS PLATE, L - 70MM - RIGHT

Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 3.5mm and locking screws.
- Plates are available in two lengths-60mm & 70mm

Indication for use

- The locking calcaneal plates address complex fractures of the calcaneus.

The locking calcaneal plate is indicated for fractures and osteotomies of the calcaneus including, but not limited to, extraarticular, intra-articular, joint depression, tongue type, and severely comminuted fractures.