



CAT. NO. SS	CAT. NO. TT	LCP COMPRESSION CANCELLOUS SCREW 4.0mm
210.120	610.520	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 30.0MM
210.125	610.525	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 35.0MM
210.130	610.530	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 40.0MM
210.135	610.535	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 45.0MM
210.140	610.540	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 50.0MM
210.145	610.545	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 55.0MM
210.150	610.550	LCP COMPRESSION CANCELLOUS SCREW 4.0mm x P.TH 60.0MM

Specification

- Available in 30mm to 60mm length with 4.0mm diameter
- Partially threaded shaft
- Available in both Stainless Steel & Titanium

Indication for use

- LCP Compression Cancellous screws are indicated for securing plates into metaphyseal and epiphyseal areas of bones
- Deep cuts and widely spaced threads enhance the holding power of the screw in soft bones

