



CAT. NO. SS	CAT. NO.TT	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM
331.003L	631.003L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 3 HOLES LEFT.
331.003R	631.003R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 3 HOLES RIGHT.
331.004L	631.004L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 4 HOLES LEFT.
331.004R	631.004R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 4 HOLES RIGHT.
331.005L	631.005L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 5 HOLES LEFT.
331.005R	631.005R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 5 HOLES RIGHT.
331.006L	631.006L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 6 HOLES LEFT.
331.006R	631.006R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 6 HOLES RIGHT.
331.007L	631.007L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 7 HOLES LEFT.
331.007R	631.007R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 7 HOLES RIGHT.
331.009L	631.009L	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 9HOLES LEFT.
331.009R	631.009R	LCP DISTAL FIBULA PLATE 2.7 MM / 3.5 MM 9HOLES RIGHT.

## Specification

- Plates(Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts LCP 2.7mm, 3.5mm cortical and locking screws.
- Plate available with 3,4,5,6,7 & 9 holes

## Indication for use

- The Wise Lock Distal Fibula Plates are intended for fixation of fractures, osteotomies and non-unions of the metaphyseal and diaphyseal region of the distal fibula, especially in osteopenic bone.