



| CAT. NO.<br>SS | CAT. NO.<br>TT | LCP DISTAL RADIUS PLATE - VOLAR PLATES                           |
|----------------|----------------|--|
| 350.503L       | 644.503L       | LCP DISTAL RADIUS PLATE - VOLAR, 3HOLES SHAFT, 5HOLES HEAD LEFT  |
| 350.503R       | 644.503R       | LCP DISTAL RADIUS PLATE - VOLAR, 3HOLES SHAFT, 5HOLES HEAD RIGHT |
| 350.505L       | 644.505L       | LCP DISTAL RADIUS PLATE - VOLAR, 5HOLES SHAFT, 5HOLES HEAD LEFT  |
| 350.505R       | 644.505R       | LCP DISTAL RADIUS PLATE - VOLAR, 5HOLES SHAFT, 5HOLES HEAD RIGHT |

## Specification

- Volar Plates are precontoured for anatomic fit on the volar aspect of the distal radius.
- Plates available in Stainless steel & Titanium material.
- Plate features locking compression holes and round locking holes that accepts 2.4mm cortical and locking screws.
- Plate available with 5 head holes & 3,5 shaft holes
- Plates are left & right configured.

## Indication for use

- It is used For the fixation of complex intra- and extra- articular fractures and osteotomies of the distal radius and other small bones.

