



CAT. NO. SS	CAT. NO. TT	LCP Distal Radius Buttress Volar Plate
350.603L	644.603L	LCP Distal Radius Buttress Volar Plate Ø 2.4mm Shaft 3 Head 5 Left
350.603R	644.603R	LCP Distal Radius Buttress Volar Plate Ø 2.4mm Shaft 3 Head 5 Right

## Specification

- Volar Plates are precontoured for anatomic fit on the volar aspect of the distal radius.
- Plates available in Stainless steel & Titanium material.
- Plate features locking compression holes and round locking holes that accepts 2.4mm cortical and locking screws.
- Plate available with 5 head holes & 3,5 shaft holes
- Plates are left & right configured.

## Indication for use

- It is used for the fixation of complex intra- and extra- articular fractures and osteotomies of the distal radius and other small bones.

