




CAT. NO. SS	CAT. NO. TT	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE
306.104L	606.104L	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 4 HOLE LEFT L 122MM
306.104R	606.104R	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 4 HOLE RIGHT L 122MM
306.106L	606.106L	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 6 HOLE LEFT L 158MM
306.106R	606.106R	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 6 HOLE RIGHT L 158MM
306.108L	606.108L	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 8 HOLE LEFT L 194MM
306.108R	606.108R	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 8 HOLE RIGHT L 194MM
306.110L	606.110L	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 10 HOLE LEFT L 230MM
306.110R	606.110R	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 10 HOLE RIGHT L 230MM
306.112L	606.112L	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 12 HOLE LEFT L 266MM
306.112R	606.112R	LCP EXTRA ARTICULAR DISTAL HUMERUS PLATE. 12 HOLE RIGHT L 266MM

Specification

- Plates (Right & Left) available in both-Titanium & Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plate available with 4,6,8,10,12 & 14 holes

Indication for use

These plates are indicated for :

- Extra-articular fractures of the distal humerus
- Malunions of the distal humerus
- Non-unions of the distal humerus

