



CAT. NO. SS	CAT. NO. TT	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX)
343.014L	643.844L	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 4 HOLES - LEFT
343.014R	643.844R	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 4 HOLES - RIGHT
343.016L	643.846L	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 6 HOLES - LEFT
343.016R	643.846R	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 6 HOLES - RIGHT
343.017L	643.847L	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 7 HOLES - LEFT
343.017R	643.847R	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 7 HOLES - RIGHT
343.018L	643.848L	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 8 HOLES - LEFT
343.018R	643.848R	LCP MEDIAL DISTAL FEMUR PLATE (TOMOFIX) 8 HOLES - RIGHT

Specification

- Plates available in Stainless steel & Titanium material.
- Plate accepts 4.5/5.0mm cortical and Locking screws.
- Plate available with 4,6,7 & 8 holes.

Indication for use

These plates are indicated for:

- Unicompartmental lateral gonarthrosis with valgus malalignment of the distal femur.
- Idiopathic or posttraumatic valgus deformity of the distal femur.