



CAT. NO. SS	CAT. NO. TT	LCP Olecranon Plate 3.5mm
349.102L	649.102L	LCP Olecranon Plate 3.5mm, 2 Hole- Left
349.102R	649.102R	LCP Olecranon Plate 3.5mm, 2 Hole- Right
349.104L	649.104L	LCP Olecranon Plate 3.5mm, 4 Hole- Left
349.104R	649.104R	LCP Olecranon Plate 3.5mm, 4 Hole- Right
349.106L	649.106L	LCP Olecranon Plate 3.5mm, 6 Hole- Left
349.106R	649.106R	LCP Olecranon Plate 3.5mm, 6 Hole- Right
349.108L	649.108L	LCP Olecranon Plate 3.5mm, 8 Hole- Left
349.108R	649.108R	LCP Olecranon Plate 3.5mm, 8 Hole- Right
349.110L	649.110L	LCP Olecranon Plate 3.5mm, 10 Hole- Left
349.110R	649.110R	LCP Olecranon Plate 3.5mm, 10 Hole- Right
349.112L	649.112L	LCP Olecranon Plate 3.5mm, 12 Hole- Left
349.112R	649.112R	LCP Olecranon Plate 3.5mm, 12 Hole- Right

Specification

- Plates available in Stainless steel & Titanium material.
- Plate accepts 3.5mm cortical and Locking screws.
- Plate available with 2,4,6,8,10 & 12 holes.

Indication for use

- LCP Olecranon Plates are indicated for fixation of fractures, osteotomies and nonunions of the olecranon, particularly in osteopenic bone.