

LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - AO TYPE STRAIGHT



CAT. NO. LVM	LCP Pelvic Recon Plate - Low Profile Ø 3.5 mm - AO TYPE STRAIGHT
305.122	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 2HOLE
305.123	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 3HOLE
305.124	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 4HOLE
305.125	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 5HOLE
305.126	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 6HOLE
305.127	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 7HOLE
305.128	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 8HOLE
305.129	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 9HOLE
305.130	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 10HOLE
305.132	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 12HOLE
305.134	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 14HOLE
305.136	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 16HOLE
305.138	LCP PELVIC RECON PLATE - LOW PROFILE Ø 3.5 MM - 18HOLE

Specification

- Available in LVM.
- Plate features locking compression holes and round locking holes that accepts 3.5mm cortical and locking screws.
- Plate available with 2,3,4,5,6,7,8,9,10,12,14,16 & 18 holes.

Indication for use

- It is used to fix Pelvic, clavicle and calcaneal fractures.