




Cat. No. SS	Cat. No TT	LCP LONG PHILOS PLATE
306.025	606.025	LCP LONG PHILOS PLATE 5 HOLE
306.026	606.026	LCP LONG PHILOS PLATE 6 HOLE
306.028	606.028	LCP LONG PHILOS PLATE 8 HOLE
306.030	606.030	LCP LONG PHILOS PLATE 10 HOLE
306.032	606.032	LCP LONG PHILOS PLATE 12 HOLE

Specification

- Plate features locking compression holes and round locking holes that accept 3.5mm cortical and locking screws
- Available in Stainless steel and Titanium

Indication for use

- Dislocated two-, three-, and four-fragment fractures of the proximal humerus, including fractures involving osteopenic bone
- Pseudarthroses in the proximal humerus
- Osteotomies in the proximal humerus
- For fractures extending into the shaft or fractures without medial support

