



Cat. No. SS	Cat. No. TT	LCP SHORT PHILOS PLATE
306.113	606.113	LCP SHORT PHILOS PLATE 3 - HOLE
306.1 14	606.1 14	LCP SHORT PHILOS PLATE 4 - HOLE
306.115	606.115	LCP SHORT PHILOS PLATE 5 HOLE

Specification

- Plate features locking compression holes and round locking holes that accept 3.5mm cortical and locking screws
- Available in Stainless steel and Titanium

Indication for use

- Dislocated two-, three-, and four-fragment fractures of the proximal humerus, including fractures involving osteopenic bone
- Pseudarthroses in the proximal humerus
- Osteotomies in the proximal humerus