

LCP POSTERIOR MEDIAL PROXIMAL TIBIA PLATE 3.5MM



CAT. NO. SS	CAT. NO. TT	LCP POSTERIOR MEDIAL PROXIMAL TIBIA PLATE 3.5MM
309.101	609.101	LCP Posteromedial Proximal Tibia Plate 3.5 mm 01 Hole
309.102	609.102	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 2 Hole
309.103	609.103	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 3 Hole
309.104	609.104	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 4 Hole
309.105	609.105	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 5 Hole
309.106	609.106	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 6 Hole
309.107	609.107	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 7 Hole
309.108	609.108	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 8 Hole
309.109	609.109	LCP Posteromedial Proximal Tibia Plate 3.5 mm 0 9 Hole
309.110	609.110	LCP Posteromedial Proximal Tibia Plate 3.5 mm 10 Hole

Specification

- Plates available in Stainless steel & Titanium material.
- Plate accepts 3.5mm cortical and locking screws.
- Plate available with 1-10 Holes.

Indication for use

- LCP Posteromedial Proximal Tibia Plate intended to buttress metaphyseal fractures of the medial tibial plateau, split-type fractures of the medial tibial plateau, medial split fractures with associated depressions and split or depression fractures of the medial tibial plateau.