



| CAT. NO.<br>SS |   | LCP PROXIMAL FEMUR HOOK PLATE         |
|----------------|---|---------------------------------------|
| 323.106        | - | LCP PROXIMAL FEMUR HOOK PLATE 6 HOLE  |
| 323.108        | - | LCP PROXIMAL FEMUR HOOK PLATE 8 HOLE  |
| 323.110        | - | LCP PROXIMAL FEMUR HOOK PLATE 10 HOLE |
| 323.112        | - | LCP PROXIMAL FEMUR HOOK PLATE 12 HOLE |
| 323.114        | - | LCP PROXIMAL FEMUR HOOK PLATE 14 HOLE |

## Specification

- Plates available in Stainless steel.
- Plate features locking compression holes and round locking holes that accepts 4.5/5.0mm cortical and Wise-lock screws.
- Plate available with 6,8,10,12 & 14 holes.

## Indication for use

**It is used to fix the following fractures of tibia :**

- Proximal shaft fractures
- Metaphyseal fractures
- Intra-articular fractures
- Periprosthetic fractures